

Consumables for Injectors and Contrast Delivery Systems



MR Suite



CT Suite



FastLoad™ CTA Dual Syringe Pack with Spikes

SKU 017355
Quantity 20/CS
Description of Product Use Syringe tray includes: (2) Disposable 200mL syringes w/ administration set, 60" low-pressure extension set, bifurcated "Y" one-way check valve and (2) spikes. For use w/ EmpowerCTA®+ Injection System.



FastLoad™ MR Syringe Pack

SKU 017356
Quantity 50/CS
Description of Product Use Syringe tray Includes: (2) 100mL syringes (1) MR connecting tube w/ bifurcated "Y" one-way check valve, (2) vented spikes for contrast & saline. For use w/ EmpowerMR® Injector System.



EmpowerCTA®+ "Y" Connecting Tubing

SKU 100114
Quantity 100/CS
Description of Product Use Administration Set for use w/ EmpowerCTA®+ Injection System. Useful when utilizing two single syringes on procedures requiring saline delivery.



FastLoad™ CTA Dual Syringe Pack

SKU 017354
Quantity 20/CS
Description of Product Use Syringe tray includes: (2) Disposable 200mL syringes w/ administration set, 60" low-pressure extension set, bifurcated "Y" one-way check valve and a fill tube. For use w/ EmpowerCTA®+ Injection System.



Imaging Bulk Pack Transfer Set

SKU 017349
Quantity 100/CS
Description of Product Use For use with Isovue® Imaging Bulk Package to transfer contrast into empty syringes in the CT Suite.



CT Exprès® Bottle Spike

SKU 640458
Quantity 60/CS
Description of Product Use For use with CT Exprès® 3D Contrast Media Delivery System. Each box contains 60 individual bottle spike bilster packs used to spike one bottle for up to 10 hours.



FastLoad™ CT Syringe Pack with Spike

SKU 017345
Quantity 50/CS
Description of Product Use Disposable 200mL syringe w/ 60" low-pressure extension set and 1 spike. For use w/ EmpowerCTA®+ Injection System.



Fluid Transfer Set

SKU 100115
Quantity 100/CS
Description of Product Use Used for transferring contrast and saline from bottles/bags to Empower branded syringes. Includes (100) 20" single use sets.



CT Exprès® Patient Set

SKU 640457
Quantity 25/CS
Description of Product Use For use with CT Exprès® 3D Contrast Media Delivery System. Each box contains 25 single-use patient sets.



FastLoad™ CT Syringe Pack

SKU 017344
Quantity 50/CS
Description of Product Use Disposable 200mL syringe w/ 60" low-pressure extension set and a fill tube. For use w/EmpowerCTA®+ Injection System.



60" Coiled Connector Tube

SKU 100113
Quantity 100/CS
Description of Product Use For use in procedures requiring an extension tubing set. Includes (100) 60" low-pressure single use sets.



CT Exprès® Day Set III HP

SKU 640460
Quantity 15/CS
Description of Product Use For use with Imaging Bulk Package and CT Exprès® 3D Contrast Media Delivery System. Each box contains 15 sets for use up to 10 hours or 4L of contrast.

To learn more about how you can place an order, how we can help you, or how you may obtain additional information about our products, please select one of the following:

Telephone: 1-877-BRACCO 9 (1-877-272-2269)

Fax: 1-866-272-1619

Email: BRACCOOTC@diag.bracco.com

Hours of Operation: M-F 8:30 AM - 6:00 PM (EST), excluding holidays



MR Suite



CT Suite

EmpowerCTA[®]+ Injector System

INDICATIONS:

The EmpowerCTA⁺ Injector System is indicated for the vascular administration of contrast and flushing media in conjunction with computed tomography (CT) scanning of the body with an optional interface to a CT scanner and an optional calculator for glomerular filtration rate (GFR).

The Extravasation Detection Accessory (EDA[™]) is an optional accessory and is indicated for the detection of extravasations of ionic and nonionic contrast during CT procedures using a power injector.

Protocols have been independently developed and are not intended as medical advice. Bracco and Bracco Injengineering S.A. shall not be responsible for any physicians' reliance on these or any other products.

CONTRAINDICATIONS:

The EmpowerCTA⁺ Injector System is not intended for use as a long-term infusion pump, nor is it intended to be used to inject any agents other than contrast or flushing media. Do not attempt to use the Injector for any other purpose (such as chemotherapy or drug infusion). The EmpowerCTA⁺ Injector System should not be used to inject substances into nonvascular body cavities. Any applications of the EmpowerCTA⁺ Injector System other than those described in its Operator Manual are inappropriate and should not be attempted.

Not all products available in all global markets.

EmpowerCTA⁺ is distributed by Bracco Diagnostics Inc.; 259 Prospect Plains Road, Bldg. H; Monroe Township, NJ 08831 USA; Phone: (800) 631-5245; Fax: (609) 514-2424; Customer Service: 1-877-BRACCO 9 (1-877-272-2269); Professional Services: 1-800-257-5181 (Option 2); EmpowerCTA⁺ Email: empowerinfo@diag.bracco.com; Website: <https://www.braccoimaging.com/us-en/products/injectors-informatics/empowercta>

EmpowerCTA⁺ is manufactured for Bracco Injengineering S.A., Avenue de Sevelin 46, CH – 1004 Lausanne, Switzerland. <http://imaging.bracco.com>

EmpowerMR[®] Injector System

INDICATIONS:

The EmpowerMR Injector System is indicated for vascular administration of contrast media and flushing media in conjunction with magnetic resonance imaging (MRI).

Protocols have been independently developed and are not intended as medical advice and Bracco shall not be responsible for any physicians' reliance on these or any other protocols.

Bracco Injengineering S.A. reserves the right at any time and without notice, to change the specifications and features described herein, or to change the production or adjust the product described.

Not all products available in all global markets.

EmpowerMR is distributed by Bracco Diagnostics Inc.; 259 Prospect Plains Road, Bldg. H; Monroe Township, NJ 08831 USA; Phone: (800) 631-5245; Fax: (609) 514-2424; Customer Service: 1-877-BRACCO 9 (1-877-272-2269); Professional Services: 1-800-257-5181 (Option 2); EmpowerMR Email: empowerinfo@diag.bracco.com; Website: <https://www.braccoimaging.com/us-en/products/injectors-informatics/empowermr>

EmpowerMR is manufactured by Bracco Injengineering S.A., Avenue de Sevelin 46 CH – 1004 Lausanne, Switzerland. <http://imaging.bracco.com>

CT Exprès[®] 3D Contrast Media Delivery System

INDICATIONS:

The CT Exprès 3D Contrast Media Delivery System is indicated for controlled automatic administration, on the venous side, of contrast media and saline, to human subjects while undergoing examination by means of a computed tomography (CT) scanner.

The CT Exprès 3D Contrast Media Delivery System is specifically indicated for use in CT procedures for the delivery of ISOVUE[®] (Iopamidol Injection) contrast media as supplied in an Imaging Bulk Package (IBP), for a maximum of 20 bottles of contrast media or a maximum of ten (10) hours, whichever comes first, per Day Set III HP disposable. The Bottle Spike disposable is for single-bottle use only and must be discarded with the contrast media bottle. The Patient Set must be discarded after each patient procedure.

The CT Exprès 3D Contrast Media Delivery System is to be used only by and under quasi-continuous supervision of trained health care professionals in an appropriate licensed health care facility, in a room designated for radiological procedures that involve intravascular administration of a contrast agent.

CONTRAINDICATIONS:

This device is not intended for injection of contrast media for coronary arteriography, or for any other use for which the device is not indicated.

Bracco Injengineering S.A. reserves the right at any time and without notice, to change the specifications and features described herein, or to change the production or adjust the product described.

Not all products are available in all global markets.

CT Exprès is distributed by Bracco Diagnostics Inc.; 259 Prospect Plains Road, Bldg. H; Monroe Township, NJ 08831 USA; Phone: (800) 631-5245; Fax: (609) 514-2424; Customer Service: 1-877-BRACCO 9 (1-877-272-2269); CT Exprès Email: empowerinfo@diag.bracco.com; Website: <https://www.braccoimaging.com/us-en/products/injectors-informatics/ct-expr-s>

CT Exprès is manufactured for Bracco Injengineering S.A., Avenue de Sévelin 46, CH – 1004 Lausanne, Switzerland. <http://imaging.bracco.com>

About Bracco Injengineering S.A.

Bracco Injengineering S.A., located in Lausanne, Switzerland, is the most recently added Business Unit of Bracco Imaging. Bracco Injengineering is committed to develop the best-in-class integrated injection solutions with a strong heritage in research and innovation. It provides quality solutions for state-of-the-art radiology centers, offering proven injection technology, built on Bracco Imaging's expertise.

Thanks to this strengthened product portfolio, including the CT Exprès 3D Contrast Media Delivery System, EmpowerCTA[®]+ Injector System, EmpowerMR[®] Injector System, NEXO[®] Contrast Management System and NEXO [DOSE][®] Multi-Modality Radiation Informatics, Bracco Imaging will be able to focus on constant innovation, not only for devices, but also for software development and data management.

CT Exprès, EmpowerCTA⁺, EmpowerMR, NEXO and NEXO [DOSE] are registered trademarks of Bracco Injengineering S.A.

EDA and FastLoad are trademarks of Bracco Injengineering S.A.

Go to: imaging.bracco.com/us-en

Bracco Diagnostics Inc. 259 Prospect Plains Road, Building H
Monroe Township, NJ 08831 USA

Phone: 609-514-2200

Toll Free: 1-877-272-2269 (U.S. only)

Fax: 609-514-2446

© 2022 Bracco Diagnostics Inc. All Rights Reserved.

Committed to Science,
Committed to You.™

